

# NEW! BRUSH BETWEEN

## WITH EASY REFILL HEADS INTERDENTAL BRUSH



**XS**  
0,4  
mm

**S**  
0,5  
mm

**M**  
0,6  
mm

**L**  
0,7  
mm

# Jordan\*

 [jordan.dentalcare.uae](https://www.instagram.com/jordan.dentalcare.uae)

 [jordan.dentalcare.uae](https://www.facebook.com/jordan.dentalcare.uae)

# Jordan\*



**XS**

**0.4 mm**



**S**

**0.5 mm**



**M**

**0.6 mm**



**L**

**0.7 mm**

PRODUCT	<b>BRUSH BETWEEN INTERDENTAL BRUSH</b>	
<b>Properties</b>	<p>The new Brush Between is an interdental brush where you can replace the brush head and keep the handle, just with a simple click! The ergonomic angled handle combined with the high-quality brush ensures superior reach and efficient cleaning. Removes plaque and food particles. The brush ensures deep cleaning to reduce and prevent inflamed gums.</p>	
<b>Find the right size</b>	<p>Available in 4 different sizes, so you can find the right fit for the spaces between your teeth. Do not force brushes between your teeth. It is also common to have multiple sizes available, as we often have different sized spaces between different teeth. If you aren't sure of which interdental brush suits your teeth best, you can get guidance from your dental hygienists or dentist.</p>	
<b>How to use</b>	<p>Insert the brush carefully into the spaces between teeth and use a gentle back and forth motion.</p>	
<b>Recommended for</b>	<p>If you have big enough spaces between teeth, braces, bridges, crowns, implants or gum diseases such as gingivitis or periodontitis.</p>	
<b>How to use refill solution</b>	<p>Press-off the used head by pressing the top of the handle together. Click-on the new head while placed in the travel case by inserting the handle, as shown in the photo.</p>	
<b>How often should the brush head be changed</b>	<p>Just as with toothbrushes, remember to switch your interdental brush regularly. Replace the brush when it starts to look worn out, to ensure continuous optimal effect and to keep your gums healthy and clean.</p>	

